

An innovative, stylish and practical addition to the new generation of EMS equipment. The Megapulse Senior 265 offers pulsed and continuous shortwave therapy with many additional features, from one of the world's largest specialised manufacturers of shortwave therapy products.

- 1,000W pulsed and 400W continuous maximum power output
- Continuous and pulsed shortwave from one unit
- Back-lit display with clear operational information
- On-screen simplicity for quick set-up
- 10 programable memory functions
- New and distinctive ergonomic design with rugged, lockable castors
- Wide range of pulse widths: 20 - 400 microseconds
- Wide range of pulse frequencies: 5 - 800 pulses per second
- Up to 30 minutes treatment time
- Choice of output modes: continuous or pulsed 1 in 3, 2 in 3 and 3 in 3
- Complete with durable electrode arms and two 100mm capacitive electrodes
- Output tester also supplied



Specification:

Mains Supply: 100-240 Vac, 1kVA max, 50/60Hz

Size: 470 x 470 x 940mm

Weight: 38kg

Frequency: 27.12MHz

Output Power: 400W maximum in Continuous Mode 1000W maximum Peak Power in Pulsed Mode

Output Modes: Continuous, Pulsed 1 in 3, 2 in 3 and 3 in 3

Classification: Class1, Type BF (EN 60601-1)

Pulse Width: 20-400 μ s

Pulse Frequency: 5-800 Pulses per Second

Timer: 0-30 Minutes with termination of output and alarm at end of treatment

Tuning: Automatic

Available Products:

Megapulse Senior 265

Megapulse Senior 265 shortwave unit including two electrode arms, two 100mm capacitive electrodes with cables, output tester and user manual

Optional Accessories:

SLA990

Monopulse inductive applicator

SLA991

Flexipulse inductive applicator

SLA993

50mm capacitive electrodes and cables-one pair

SLA996

One pair rubber electrodes with four felt spacers (180 x120mm)

SLA9960

Felt spacers for SLA9960

SLA997

One pair rubber electrodes, 260 x 180mm with four felt spacers

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk