

The NeuroTrac Simplex is an advanced and simple to use single channel EMG biofeedback unit, featuring a clear display for the patient.

- Single channel EMG biofeedback
- Structured protocols for patient progress reports (optional software)
- Measurements from 0.2?V to 2000 ?V
- Use for applications such as treatment for incontinence, shoulder and neck pains, patello femoral pain, repetitive strain injury, stroke and much more.
- Monitor patient usage between treatments
- Performance software available

“EMG Biofeedback) acts as an enhancer of the therapy, enabling the patient (and the therapist) to make more effective and rapid progress towards the rehabilitation goal”

Professor Tim Watson Electrotherapy on the web



Specification:

Size: 140 x 90 x 29mm

Weight: 0.2kg

Power: 9v Battery PP3 (6LR61)

EMG, Stimulation: TENS and ECS

Outputs: For Full Details See Operating Manual

Display: 100 x 65 Pixel LCD with Backlight

Optional Accessories:

RB495

Femelex vaginal probe-single use (26mm)

EMS546

NeuroTrac software

EMS547

30mm diameter self-adhesive electrodes

Available Products:

NeuroTrac Simplex EMG Biofeedback unit

NeuroTrac Simplex EMG Biofeedback unit

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk