

For outside EU only.

A small wax bath ideal for hand therapy. The temperature on this unit can be adjusted to suit the needs of the patient and is thermostatically controlled.



## Specification:

**Power Requirement:** 200/240V 50/60Hz, (120V 60Hz) 40/90W

**Wax Capacity:** 2.7kg

**Dimensions:** 457 x 311 x 216mm

**Weight:** 6.8kg with Wax

**Classification:** Class-1, Type B (EN 60101-1)

## Available Products:

### Small Wax bath

Para-Care wax bath (220-240V), 120V also available on request.

## Optional Accessories:

### T1A

Wax bath thermometer T1A (supplied as standard with LNC884 and LNC882)

### LNC885

Tubular steel stand for wax bath (supplied flat-packed)

### P11

Paraffin wax 25kg

## Contact Information:

### UK Contact

EMS Physio Ltd  
Grove Technology Park  
Downsview Road  
Wantage  
Oxfordshire OX12 9FE  
England

T +44 1235 772272

F +44 1235 763518

E [sales@emsphysio.co.uk](mailto:sales@emsphysio.co.uk)