

Allows up to 10kg of wax to be purified at any one time. The tapered shape of the wax container ensures easy wax removal after cooling.



Specification:

Power Requirement: 200/240 Vac. 2.3kW

Water Capacity: 10.5 litres

Wax Capacity: 10kg

Dimensions Outer Unit: - 306mm Diameter x 525mm, Inner Unit - 300-225mm Tapering x 275mm Deep

Temperature Range: 0-100 °C

Certification: EMC Directive and LVD

Available Products:

EMS Wax purification unit

Wax purification unit (220-240V) includes tapered wax container (LNC887)

Optional Accessories:

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk