

The driving force behind the MEDI-LINK system. This controls and monitors each MEDI-LINK module, displaying the details of each therapy on a bright, back-lit screen. The keypad enables each therapy to be chosen and treatments to be selected.

Once a MEDI-LINK control module is purchased, simply choose the modules to build a custom-made therapy centre.

- Clear set-up and treatment selection on a bright back-lit screen
- Multilingual display available from the system set-up menu
- Touch sensitive keypad with raised keys
- Compact and lightweight
- Flip-up screen can be locked down for transportation



Specification:

Power Input: 100-240Vac 200va

Display: 240 X 128 Pixel Graphics Backlit LCD

Classification: Class 1 BF (EN60601-1)

Size: 240 x 210 x 100mm

Weight 2 kg:

Available Products:

MEDI-LINK Control Module

MEDI-LINK control unit with carrying handle

Optional Accessories:

EMS158

Primo trolley

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk