

The MEDI-LINK dual-frequency ultrasound module has two sockets for permanent connection of both large and small treatment heads - saving time. A simple press of the key on the control module switches between treatment heads. Clear indications on the control module screen show selections of time, treatment head, mode of output and contact.

- Large or small treatment heads available
- Dual-frequency output through each of the large or small treatment heads
- Light indicator in the treatment head clearly shows successful contact
- Treatment time stops with visual indication if patient contact is insufficient
- Wide range output of modes (1:1, 1:2, 1:4, 1:9, continuous and user-defined)
- Can be connected to the MEDI-LINK interferential module for combination therapy
- Memory function for easy recall of favourite programs

"One of the therapeutic effects for which ultrasound has been used is in relation to tissue healing"

Professor Tim Watson Electrotherapy on the web



Specification:

Frequency: 1MHz And 3MHz

Treatment Time: 0-30 Minutes (Treatment Linked)

Intensity: 0-2W/cm² Modes CW And Pulsed 1:1, 1:2, 1:4, 1:9 And User-defined

Contact Monitor: Luminous Ring on Treatment Head Large Transducer ERA 5cm² BNR

<6 Small Transducer ERA 0.5cm² BNR <6

Size: 210 x 80 x 100mm

Weight: 1kg

Available Products:

MEDI-LINK ultrasound module

MEDI-LINK dual-frequency ultrasound module (must be used with MEDI-LINK control unit)

Optional Accessories:

ML9120

Large treatment head 5 sq cm

ML9130

Small treatment head 0.5 sq cm

ML9140

Combination lead

ML9150

Large angled treatment head shower style 5 sq cm

EMS502A

1 litre of gel

EMS502

8X250ml bottles of gel

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk