

The MEDI-LINK EMG biofeedback module offers patient education and rehabilitation on a wide variety of surface muscles, as well as continence re-education - from just one unit. With an on-screen graph for instant evaluation, this system is highly sophisticated yet easy to operate.

- Single or dual-channel quantitative information
- Suitable for a wide range of muscle re-education conditions - anterior knee, unstable shoulder, continence re-education, oral pharyngeal dysphasia, TMJ and bruxism
- Ideal for rehabilitation in post-operative trauma, teaching muscle relaxation, posture improvement, phantom limb pain and goal setting
- Real time signal displays in both bar chart and graph form
- Memory of previous settings - saving time

"EMG Biofeedback acts as an enhancer of the therapy, enabling the patient (and the therapist) to make more effective and rapid progress towards the rehabilitation goal"

Professor Tim Watson Electrotherapy on the web



Specification:

No. of Channels: 2

Display: Line and Bar Graphs and Numeric Value

Sensitivity: 20 μ V To 2 mV, Full Scale in Steps

Audible Response Alarm: Above or Below Threshold from MEDI-LINK Speaker

Frequency Response: 20-500Hz (Wide) or 100-200Hz (Narrow) +50/60Hz Notch

Size: 210 x 80 x 100mm

Weight: 1kg

Available Products:

MEDI-LINK EMG biofeedback module

MEDI-LINK EMG biofeedback module (must be used with MEDI-LINK control unit)

Optional Accessories:

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk