

Robust two-section couch ideally suited for continuous use in areas such as out patient departments, gp surgeries, colleges, private clinics etc.

Features:

- unique lifetime guarantee
- fail-safe backrest locking mechanism
- low running costs
- hardwearing antimicrobial vinyl
- able to lift 160kg from its lowest height
- electric or hydraulic operation

**Specification:**

Width: 65cm

Length: 189cm

Height range: from 46cm to 95 cm

Maximum lifting weight: 160kg (25 stone)

Breathing hole: option

Packed weight: 82kgs

Shipping dimensions: 190 x 66 x 64cm

Available Products:**2 Section Hydraulic Plinth RH075A**

Two-section hydraulic plinth

2 Section Hydraulic Plinth RH075B

Two-section hydraulic plinth - with breathing slot

2 Section Electric Plinth RH076A

Two-section electric plinth

2 Section Electric Plinth RH076B

Two-section electric plinth - with breathing slot

Optional Accessories:**RH088**

Paper roll holder

RH092

Gynaecology stirrups - chrome

RH093

Foot switch

RH094

Cot side (per side)

RH096

Towelling cover

Contact Information:**UK Contact**

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emspphysio.co.uk