

Very popular versatile three-section couch, firm and stable providing the ideal base for general treatment techniques and heavy manipulation. Commonly used in areas such as physiotherapy, GP surgeries, alternative therapies, private clinics etc.

Features:

- unique lifetime guarantee
- fail-safe backrest locking mechanism
- low running costs
- hardwearing antimicrobial vinyl cover
- able to lift 160kg from its lowest height
- electric or hydraulic operation.

**Specification:**

Width: 65cm

Length: 189cm

Height range: from 46cm to 95 cm

Maximum lifting weight: 160kg (25 stone),

Breathing hole: option

Packed weight: 82kgs

Shipping dimensions: 190 x 66 x 64cm

Available Products:**3 Section Plinth - hydraulic RH078A**

Three-section plinth -hydraulic

3 Section Plinth - hydraulic RH078B

Three-section plinth - hydraulic with breathing slot

3 section plinth - electric RH079A

Three-section plinth - electric

3 section plinth - electric RH079B

Three-section plinth -electric with breathing slot

Optional Accessories:**RH088**

Paper roll holder

RH092

Gynaecology stirrups - chrome

RH093

Foot switch

RH094

Cot side (per side)

Contact Information:**UK Contact**

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emspphysio.co.uk