

This neurology 'double' plinth offers all the features of our standard Wide Neurology plinth with the added benefit of being able to lift a very heavy patient to a suitable working height when the need arises. The standard platform is 102 cm although wider versions can be made to order. Motorised backrest recommended.

Features: unique lifetime guarantee, fail-safe backrest locking mechanism, low running costs, fully joist accessible, hardwearing antimicrobial vinyl, able to lift 320 kg from its lowest height, electric operation.



Specification:

Electric height adjustment

Width: 102cm

Length: 190cm

Height range: from 48cm to 99cm

Maximum lifting weight: 160kg (25 stone)

Breathing hole and plug fitted as standard

Packed weight: 130kg

Shipping dimensions: 190x102x57cm

Available Products:

Heavy-duty electric neurology plinth RH081H

Heavy-duty electric neurology plinth

Optional Accessories:

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk