

Friction free rolling top traction table available with or without RH087 traction machine.
Supplied with all accessories necessary for both lumbar and cervical traction.

Features:

- Unique lifetime guarantee
- Low-running costs
- Fully hoist accessible
- Hardwearing antimicrobial vinyl
- Able to lift 160kg from its lowest height
- Electric or hydraulic operation



Specification:

Electric or hydraulic options

Width 65cm

Height 65cm

Length 199cm

Height range 46cm to 95cm

Maximum lift 160kg (25 stone)

Breathing hole and plug fitted as standard

Available Products:

Traction Table, hydraulic

Friction free hydraulic rolling top traction table. Supplied with all accessories necessary for both lumbar and cervical traction.

Traction Table, electric

Friction free electric rolling top traction table. Supplied with all accessories necessary for both lumbar and cervical traction.

Optional Accessories:

RH086

Flexion stool

RH088

Paper roll holder

RH092

Gynaecology stirrups - chrome

RH093

Foot switch

RH094

Cot side (per side)

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emspphysio.co.uk