

An easy-to-assemble, lightweight portable plinth which is also very strong and rigid. It has a breathing slot and adjustable head rest, with a carrying bag supplied as standard. Optional variable height and extra width.



Specification:

Weight 18 kg
Width 61cm
Length 183cm
Height 76cm
Max User Weight 25 stone

Optional Accessories:

RH113
Portable plinth carrying bag

Available Products:

Portable Plinth

An easy-to-assemble, lightweight portable plinth. It has a breathing slot and adjustable head rest, with a carrying bag supplied as standard.

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272
F +44 1235 763518
E sales@emsphysio.co.uk