

The new Primo Therasonic 460 offers you an impressive range of choices for your ultrasound treatment programme:

- Dual frequency ultrasound (1MHz and 3MHz)
- Programmable memory
- Pulsed or continuous ultrasound output
- Ergonomic treatment head

What's more, Primo offers some innovative new features:

- Touch screen technology
- Easier to clean-reducing potential cross infection
- Simple to set up
- Easy to use
- Lightweight and portable
- Mains powered as standard
- Optional rechargeable battery available

"One of the therapeutic effects for which ultrasound has been used is in relation to tissue healing"

Professor Tim Watson Electrotherapy on the web



Specification:

Size (h/w/d) 90x230x290mm

Weight 1.3 kg

Treatment Time 0-30 min

Treatment Programs 10 user-defined set-ups

Ultrasound Frequency 1MHz and 3MHz

Maximum Intensity 1.5 W/cm² in CW, 3.0 W/cm²

Maximum Output Power 6W average (Mains operated only)

Output Modes CW pulsed 1:1, 1:2, 1:4 and 1:9

Pulse Duration 2ms

Contact Monitor Light on transducer or audio signal Transducer ERA 4 cm²

Available Products:

Primo Therasonic 460 Mains EMS460

Primo Therasonic 460 Dual frequency ultrasound unit mains only (1 MHz and 3MHz) complete with large treatment head, one bottle of coupling gel (170ml)

Primo Therasonic 460 Mains/ Battery EMS460B

Primo Therasonic 460 Dual frequency ultrasound unit with battery option (1 MHz and 3MHz) complete with large treatment head, one bottle of coupling gel (170ml)

Optional Accessories:

EMS158

Primo trolley

EMS530

Primo carrying case

EMS502A

1 litre of gel

EMS502

8X250ml bottles of gel

SLA9160

Combination lead

SLA9130

Small Treatment Head (for use with EMS460, EMS455 & EMS855)

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk