

A complete work station for the ultimate exercise programme, this multigym can accommodate several occupants at the same time and has a leg press and dip station - ideal for fitness and rehabilitation.

- Space-saving leg press
- Adjustable in length
- Weight load 130kg (can be increased to 160 kg)
- Dips and intensive stomach-muscle training (leg lifts)
- Adjustable leg stretcher
- Seat adjustable in height
- High-quality padded seat and back rest
- 80 kg weight block, load variable from 5-80 kg (can be increased to 100 kg)



Specification:

Maximum load: 150 kg (user weight)
Size: L = 200cm, W = 200 cm, H = 200cm
Delivery & Installation recommended

Optional Accessories:

Available Products:

Kettler Ultra Multigym

A complete work station for the ultimate exercise programme

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272
F +44 1235 763518
E sales@emsphysio.co.uk