


Performa Knee CPM including accessories Range of Motion, Extension Flexion

- Knee -3° 130°
- Hip 8° 80°
- Ankle Plantar flexion 40° Dorsiflexion 30°
- Internal rotation 30°, External rotation 30°



Specification:

Patient height: 112cm to 206cm Speed from 50°/minute to 220°/minute

Weight: 15kg

Length: 109cm

Width: 33cm

Height: 35cm

Voltage: from 100 to 240 volts, 50/60Hz

Electrical power: 60VA

Available Products:

Performa Knee CPM

Performa Knee CPM including accessories Range of Motion, Extension Flexion

Optional Accessories:

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk