

Maestra portable hand CPM including accessories

- Range of motion: 15° hyperextension to 270° flexion (composite fist)
- Speed (5 gradations): 27 to 135 seconds (cycle of 0°-270°-0°)



Specification:

Available Products:

Maestra portable hand CPM

Maestra portable hand CPM machine including accessories

Optional Accessories:

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk