

Pulse Press Multi 3 Pro, 3 chamber sequential gradient compression. Pulse Press products are manufactured in the UK, to the highest standard. The variety of pulse presses are powerful, hi-tech and have failsafe automatic deflation.



Specification:

Pressure 30-100mmHg

Compressor 14-20 ltrs/min

Cycle Times:

Inflation 5-120 sec

Deflation 5-60 sec

Garments 1ch & 3 ch

Size 166x135x150mm

Weight 1.85kg

Supplied with 3-to-1 chamber socket for solo garments

Available Products:

Pulse Press Multi 3 Pro

Pulse Press Multi 3 Pro, 3 chamber sequential gradient compression.

Optional Accessories:

PCT011

Half Arm- single chamber

PCT012

Long arm sleeve - single chamber

PCT013

Wide long arm sleeve - single chamber

PCT014

Foot and ankle legging - single chamber

PCT015

Half legging - single chamber

44-1

Long leg- single chamber

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk