

The SoLo Multidyne 965 provides all low and medium frequency outputs from one unit.

- Stimulation on the SoLo Multidyne 965 includes Interferential (2 and 4-pole), Russian, Diadynamic, TENS, Sinusoidal, Faradic, Galvanic, Interrupted Galvanic, Trabert and Medi-Wave
- Lightweight and portable with optional mains/battery function
- Battery power indicator on screen
- Three carrier frequencies - 2kHz, 4kHz and 8kHz
- Vector facility
- 0-100mA peak constant current or 0-70/140V peak constant voltage, depending upon mode
- Provides combination therapy when linked to EMS SoLo ultrasound units
- Operates in conjunction with EMS vacuum unit, if required
- Memory function for easy recall of favourite programs
- Optional shoulder bag and treatment trolley available



Specification:

Power Input: 100-240 Vac 50VA 50/60Hz

Output Voltage: 0-60V, 0-70v Peak or 0-140V Peak Dependent on Mode

Display: LCD, Size 240 x 210 x 100mm

Battery: (Optional) Internal Rechargeable (Nimh)

Weight: 1.5kg.

Output Type: Constant Current or Constant Voltage

Treatment Time: 0-30 Minutes

Output Current: 0-70mA or 0-100mA Dependent on Mode

Treatment Programmes: 10 User-defined

Individual Waveforms: See Operating Manual for Full Details

Output Mode: 2 or 4-pole

Available Products:

SoLo EMS965 Stimulation Unit

Mains only unit with 2/4-pole patient lead and our electrode connection cables (blue and yellow), four medium (100x70mm) rubber electrodes, four medium sponge covers and two stretch bandages (1200x75mm)

SoLo EMS965 Stimulation Unit

Mains/battery unit with 2/4-pole patient lead and our electrode connection cables (blue and yellow), four medium (100x70mm) rubber electrodes, four medium sponge covers and two stretch bandages (1200x75mm)

Optional Accessories:

NC3052A

Small sponge electrode covers-set of 4

NC3052B

Small conductive rubber electrodes -set of 4 (70 x 50mm)

NC3053A

Medium sponge electrode covers-set of 4

NC3053B

Medium conductive rubber electrodes-set of 4 (100 x 70mm)

NC3054A

Large sponge electrode covers-set of 4

NC3054B

Large conductive rubber electrodes-set of 4 (130 x 100mm)

Contact Information:

UK Contact

EMS Physio Ltd
Grove Technology Park
Downsview Road
Wantage
Oxfordshire OX12 9FE
England

T +44 1235 772272

F +44 1235 763518

E sales@emsphysio.co.uk