



EMS Physio Ltd

Manipulation Couch

An extremely strong and rigid 3 section couch designed for various examination and manipulation procedures and suitable for constant heavy use. This model features a twin-spine frame, providing additional strength, rigidity and long life. Either section can be raised up to 85° to fully support the patient in sitting or they can be lowered to -25° for prone treatment procedures. One end has a breathing hole and plug for this purpose.




Standard Features

- UK manufactured
- Lifetime frame & 5 year parts guarantee
- 16 year average lifespan
- Choice of electric or hydraulic operation
- Very strong & rigid twin-spine frame
- Breathing hole & plug
- Hinge covers
- Gas strut operated head & foot sections
- Chrome release levers
- Retractable wheel base
- Adjustable foot for uneven floors
- Hoist clearance 13cm (5")

Overall Dimensions

Length:	200cm (79") approx 200-225cm (79"-89") with foot extension option
Width:	63cm (25") approx 70cm (28") option available
Height:	47cm - 97cm (18" - 38") approx



Options & Accessories



Footswitch



Dual Footswitch



Foot Ring Switch



Paper Roll Holder



Foot Extension





Linked Castor Base

NHS Approved Supplier



t: +44 (0) 1235 772272
e: sales@emsphysio.co.uk
www.emsphysio.co.uk

Operational Parameters



EMS Physio Ltd order codes
COU013 Manipulation couch

For Your Piece of Mind

- Lifetime frame and 5 year parts (accidental damage to upholstery excluded)
- Optional extended warranty to 10 years
- Average life span of 16 years
- Designed and manufactured in-house in the United Kingdom
- CE marked
- Externally audited to ISO9001:2015 by LQRA

Standards

- IP Rating - IPX6
- EE Classification - EU MDR 2017/745 Annex VIII Class 1
- Electrical Standard - IEC 60601-1
- Vinyl flame retardant to CRIB 5 (BS 5852)
- Foam flame retardant to BS 5852: Part 2: 1982, Ignition Source 5 [CRIB 5]

Delivery Promise

- Delivery to UK mainland in just 7-10 working days
- Faster service maybe available if required.

Upholstery Vinyl

- Flame retardant to CRIB 5
- Patented anti-microbial technology proven to inhibit the growth of E-Coli and MRSA
- Excellent levels of wear and longevity
- Fully washable and impermeable, to aid infection control
- Suitable for cleaning and disinfecting with dilute solutions of Hypochlorite (10,000ppm)
- Please note; Unless specified Sapphire upholstery will be supplied



Jasmine



Battleship



Cobalt



Cocoa



Grape



Mulled Wine



Pillarbox



Gingersnap



Candy



Applemint



Sapphire



Atlantic Blue



Lupin



Denim



Cool Blue



Rainforest



Wasabi



Citrus Green



Marigold



Almond

EMS PHYSIO Ltd
Downsview Road
Wantage
Oxfordshire OX12 9FE

t: +44 (0) 1235 772272
f: +44 (0) 1235 763518
e: sales@emsphysio.co.uk
www.emsphysio.co.uk