

Wheely Stool

This wheely stool has a comfortable 41cm diameter seat and a 63cm 5 castor base offers maximum stability for this gas lift chair which is available in a range of colours to match your couch.

Supplied with a 1 year unconditional guarantee.

Technical Specification

Seat Diameter 410mm

Overall Height 420 – 560mm

Diameter of Base 630mm

Ordering Information

Wheely Stool order code **WS**

Upholstery colours

High quality upholstery is a must for any good treatment couch being the interface between the couch and patient it must remain clean and resistant to damage. We use only the best quality vinyl proven over many years and all feature an anti-microbial additive to help prevent cross infection.

All of our products are available in the upholstery colours shown

